Osteosynthese/ Osteosynthesis

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Instructions for use

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PTI0029-004-0322
(Ausgabedatum: 21.04.2022)

Inhaltsverzeichnis / Content

Deutsch	4
English	35
Español	63
Português	94
Italiano	125
Français	155
USA English	186
USA Español	210

Inhalt

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2.1. Intended Purpose

WARNING

Use of implants contrary to intended purpose

- Risk of injury due to implant failure!
- Implants must only be used in accordance with intended purpose.

2.1.1. Foot Plate Systems

- The plate-and-screw systems are implants used to perform and stabilize an osteotomy and/or osteosynthesis, including (among other things) to correct deformities of the forefoot, midfoot, or hind foot.

2.1.2. Screws

- Screws are implants used to secure bone fractures, stiffen joints (arthrodesis) or reconstruct bones (osteotomies) in the forefoot, midfoot, or hind foot.

2.1.3. Staples

- Staples serve to correct misalignments in the forefoot region in relation to the presence of a hallux valgus interphalangeus.

2.1.4. Kirschner Wires

- Kirschner wires facilitate the osteosynthesis process, help secure bone fractures, and guide instruments/implants in the fore-, mid-, or hind foot region.

2.2. Indications

2.2.1. Foot Plate Systems

MetaFix™ MTP

- Hallux rigidus
- Revision following unsuccessful Keller-Brandes operation
- Revision following unsuccessful arthroplasty

MetaFix™ I

- Hallux rigidus
- Hallux valgus
- Arthrosis of the tarsometatarsal (TMT) I joint
- Tarsometatarsal (TMT) I joint hypermobility

MetaFix™ BLP

- Hallux rigidus
- Hallux valgus
- Arthrosis of the tarsometatarsal (TMT) I joint
- Tarsometatarsal (TMT) I joint hypermobility

MetaFix™ MR10

- Hallux rigidus
- Hallux valgus

MetaFix™ MTPlus

- Hallux rigidus
- Revision following unsuccessful Keller-Brandes operation
- Revision following unsuccessful arthroplasty

MetaFix™ OpenWedge

- Hallux valgus

MetaFix™ Plantar BG 10

- Hallux valgus
- Arthrosis of the tarsometatarsal (TMT) I joint
- Tarsometatarsal (TMT) I joint hypermobility
- Metatarsal I elevatus deformities

MetaCun™

- Arthrosis of the tarsometatarsal (TMT) joints (Lisfranc joints)
- Arthrosis of the cuneonavicular joint
- Arthrosis of the calcaneocuboid joint

DuoMetaCun™

- Arthrosis of the tarsometatarsal (TMT) joints (Lisfranc joints)
- Arthrosis of the cuneonavicular joint
- Arthrosis of the calcaneocuboid joint
- Arthrosis of the talonavicular (TN) joint

MetaStep™ Calcaneus

- Acquired adult flatfoot deformity (AAFD) or tibialis posterior tendinitis
- Pes cavus deformities
- Rigid pes planus

MetaFix™ TTF

- Late-stage arthritis of the upper subtalar joint (including rheumatoid arthritis, osteoarthritis, post-traumatic arthritis)
- Arthritic varus misalignment
- Arthritic valgus misalignment
- Revision following unsuccessful subtalar joint endoprosthesis
- Revision following unsuccessful subtalar joint arthrodesis

2.2.2. Screws

DuoThread™ Scarf N

- Hallux Valgus
- Hallux Rigidus
- Arthrosis of the metatarsophalangeal joint (MTP joint)
- Interphalangeal arthrosis

Merete® Cannulated HCS

- Hallux Valgus
- Hallux Rigidus
- Arthrosis of the metatarsophalangeal joint (MTP joint)
- Fixation of intra-articular and extra-articular fractures and non-unions of small bones and bone fragments (incl. calcaneus, talus, metatarsus) in the mid- and hind foot region

TwistCut™ Snap-Off

- Distal oblique midfoot osteotomy (Weil osteotomy)
- Distal chevron osteotomy
- Monocortical fixation of small bone fragments

Mecron™ Screw Sytem

- Bone reconstructions
- Osteotomies
- Arthrodeses
- Joint fusions
- Band fixations
- Fracture repairs and fixations

2.2.3. Staples

- Hallux valgus interphalangeus

2.2.4. Kirschner Wires

- Bone fragment fixation and guiding of instruments/implants
- Various osteosynthetic procedures, fixation and fracture treatments involving small fragments

2.3. Contraindications

- Acute or chronic infections of the foot or immediate vicinity
- Allergies to any of the materials used

2.4. Factors Interfering with success

- Obesity and pre-obesity
- Local bone tumours
- Osteoporosis or osteomalacia
- Systemic or metabolic disorders
- Alcohol or drug abuse
- Physical activities associated with strong shocks which could result in the implant being exposed to impacts and/or excessive loads (e.g. hard physical labour, certain types of sport)
- Patients with pathological mental or neurological conditions, or patients who are not capable of following the necessary post-operative care instructions

2.5. Possible Side Effects

WARNING

Examination of patient using MRI

- Risk of injury due to alternating magnetic fields!
- Merete GmbH does not authorise the use of MRI examinations in conjunction with the components described in these user instructions.
- Always perform an individual risk-benefit analysis.
- Check whether other imaging procedures can be used to achieve the desired diagnostic goal.

The side effects listed below are among the most common adverse effects of implantation procedures:

- Loosening and dislocation of the implant due to unfavourable anchoring/suboptimal positioning in the bone
- Implant failure due to excessive stress
- Early and late onset infections
- Wound haematomas and delayed wound healing
- Vascular disorders, including vein thrombosis, pulmonary embolism and heart failure
- Heterotopic ossification
- Nerve damage as a result of surgical trauma
- Difficult removal of the implant with long period of implantation (> 6 months) due to grown in bone
- Pain (chronic regional pain syndrome, hardware-related)
- Pseudarthrosis / Insufficient healing of bone fracture sites
- Delayed fracture healing
- Misaligned fracture healing
- Improper fixation
- Implant failure, implant breakage and/or metal wear
- Loss of correction
- Soft tissue irritation
- Pressure
- Scar infection and pain
- Incision dehiscence
- Neuritis
- Under- or overcorrection
- Uncontrolled fracture of the lateral corticalis
- Shortening of the first metatarsals
- High recurrence rate of hallux valgus in growing patients

- Transfer metatarsalgia
- Deformity of the metatarsal I
- Residual hallux interphalangeus
- Rupture of the tibialis anterior tendon
- Necrosis (avascular and soft tissue)
- Extremity Amputation
- Ulcer
- Skin Numbness and inflammation

Forefoot implants

- Impaired function of the metatarsophalangeal I joint
- Osteonecrosis of the first metatarsal head
- Arthrosis of the metatarsophalangeal I joint
- Lesion of the extensor or flexor tendon resulting in functional impairment
- Painful osteophyma of the metatarsophalangeal joint
- Contracture of the interphalangeal joint

Midfoot implants

- Tendinitis of the M. tibialis anterior
- Tendinitis of the M. extensor hallucis longus
- Capsulitis of the metatarsophalangeal II joint
- Sesamoidosis

Hindfoot implants

- Calcaneodynia
- Peroneal impingement syndrome
- Traumatic arthritis

2.6. System Compatibility

Implants may only be combined with other Merete GmbH products. Implants may only be combined with product-specific system components as per Merete specifications.

Implant Materials

Materials used in manufacturing implants are indicated on the label:

Foot Plate Systems

Foot plate systems are made of the following material:

- TiAl6V4 ELI (DIN EN ISO 5832-3)

Screws

Screws are made of the following material:

- TiAl6V4 ELI (DIN EN ISO 5832-3)

Staples

Staples are made of the following materials:

- TiAl6V4 ELI (DIN EN ISO 5832-3)
- X2CrNiMo 18-15-3 (DIN EN ISO 5832-1)

Kirschner Wires

Kirschner wires are made of the following materials:

- TiAl6V4 ELI (DIN EN ISO 5832-3)
- X2CrNiMo 18-15-3 (DIN EN ISO 5832-1)

Additional information on the chemical composition and mechanical properties of the materials used is available from Merete GmbH on request.

3. Use

3.1. General Instructions

WARNING

Use of damaged or defective Instruments/Implants

- Risk of injury due to premature implant failure!
- Instruments/Implants with identifiable damage may not be used.
- Avoid notches, scratches or bending of the implant in order to preserve its stability.

WARNING

Use of implants contrary to intended purpose

- Risk of injury due to implant failure!
- Implants must only be used in accordance with intended purpose.

WARNING

Combination with products from other manufacturers

- Risk of injury due to implant failure!
- Do not combine implant components with products from other manufacturers.

WARNING

Foreign bodies (e.g., tissue, bones) between implant components

- Risk of injury due to implant failure!
- Thoroughly clean any foreign bodies from implant components.

WARNING

Use of Implant/Instrument contrary to intended purpose

- Damage to/destruction of instrument/implant and injury to patient!
- Ensure correct implant/instrument handling. Do not misuse!

Prior to using Merete products, surgeons and assisting staff must study in detail the safety information in this product information sheet as well as the product-specific guidelines listed in the surgical technique. The relevant documentation is available from Merete on request. Surgeons must also be aware of any remaining risks associated with the products he or she intends to use, and must inform patients of such risks in advance. Implant operations must only be performed by surgeons who are not only qualified to carry out such operations, but also have extensive knowledge and experience in this field. The surgeon bears all responsibility for adverse effects or complications arising from misdiagnosis, improper surgical technique, incorrect implant selection or handling, or failure to observe the safety instructions provided in this product information sheet. Neither the manufacturer nor authorised Merete product representatives may be held liable in such cases. Before operating, study the techniques outlined in the surgical technique carefully. Training in the correct handling of implants and instruments by an authorised Merete representative is essential.

Patients must be informed about any alternative treatment methods and adequate systems having the same indication. The appropriate Merete instruments must be used for both implantation and explantation. The use of Merete instruments for purposes other than those intended is not permitted. Explanted components are not to be reused under any circumstances. Products which are supplied sterile may not be resterilised. These products are not known to interact with any medications.

For further information, consult the catalogues and operations instructions for the products in question. In case of doubt, contact your sales representative or the Merete GmbH Customer Care/Sales team prior to operating.

3.2. Instructions on Patient Information

Surgeons must make their patients aware of information mentioned in Section 2, i.e. any factors that could hinder the success of the operation as well as potential complications that may occur as a result of specific indications. There is a risk of the implant failing prematurely, or loosening due to insufficient osseointegration, changes in load bearing behaviour, and/or tissue reactions to the implant and its abrasion particles; see also 2.5. Further risks exist in connection with the general surgical procedures to be performed, as well as risks of an operation under anaesthesia. Patient has to be informed about alternative methods of treatment and other systems with adequate intended uses.

Patients also need to be informed of any measures they can take to help prevent potential consequences of these factors. The operating surgeon should keep a written record of all information given to the patient.

3.3 Features

WARNING

Combination with products from other manufacturers

- Risk of injury due to implant failure!
- Do not combine implant components with products from other manufacturers.

WARNING

Combination of implant components of different sizes

- Damage to implant components!
- Combine only components of the same size.

WARNING

Use of instruments with electrical energy

- Risk of injury due to implant failure!
- Do not damage the surfaces of the implants under any circumstances.

The osteosynthetic implants have not been evaluated in terms of safety and tolerability within MR environments. They have not been tested for warming, migration or image artefacts in MR environments. The safety of the osteosynthetic implants in MR environments is not known. MR examinations on patients with the implants could result in injuries.

3.3.1. Foot Plate Systems

Fixed-angle foot plate systems consist of implants (plates and screws) and instruments (trial plates), and are used in the osteosynthesis of small fragments. The systems' plates and screws are made of TiAl6V4 ELI (DIN EN ISO 5832-3), and are available in a variety of sizes and lengths. They may only be used in combination with other Merete products belonging to the same system.

3.3.2. Screws

The screws are used to fix small bone fragments and stiffen small joints. All osteosynthetic screws are made of TiAl6V4 ELI (DIN EN ISO 5832-3) and are available in a variety of sizes and lengths. They may only be used in combination with other Merete products belonging to the same system.

3.3.3. Kirschner Wires

Kirschner wires are used to secure bones, bone fragments, instruments, or implants to bones, either temporarily or permanently. The wires are practically universal in their range of applications, and can be used with a variety of instruments and implants. K-wires are made of X2CrNiMo (DIN EN ISO 5832-1) or TiAl6V4 ELI (DIN EN ISO 5832-3).

3.3.4. Staples

Staples serve to correct misalignments in the forefoot region in relation to the presence of a hallux valgus interphalangeus. The staples are made of TiAl6V4 (DIN EN ISO 5832-3) or X2CrNiMo (DIN EN ISO 5832-1).

3.4. Pre-Operative Planning

WARNING

Combination of implant components of different sizes

- Damage to implant components!
- Combine only components of the same size.

WARNING

Combination with products from other manufacturers

- Risk of injury due to implant failure!
- Do not combine implant components with products from other manufacturers.

Note: The following warning applies only to the foot plate systems.

WARNING

Implantation of trial implants

- Risk of injury due to breakage of trial implant!
- Only use trial implants in order to select a suitable permanent implant.
- Trial implants are not suitable for permanent implantation.

Operation planning should be done on the basis of in-depth evaluations of patient X-rays, which provide the information necessary for determining the appropriate prosthesis type, size and possible combinations. The size of the osteosynthetic implant is selected based on the given anatomical conditions. It is also important to carry out pre-operative tests on the patient to rule out allergic reactions to the implant materials. Use trial components during the operation in order to check implant size and positioning. The components planned for implantation must be on hand in all available sizes. It is absolutely essential to check the implant label (type, item number, material and size) against the information on the package before insertion. Use the additional enclosed labels for documentation in surgery reports and for reordering. Check that all implant surfaces are free of soiling.

Important: Always check the sterilisation expiration date!

When combining with other Merete components, observe the advice in the relevant user instructions, instrumentation guides and surgery videos, if any. In case of doubt, contact Merete GmbH.

3.5. Implant Handling

3.5.1. Instructions for Handling

WARNING

Use of implants which have been previously used

- Risk of injury due to premature implant failure!
- Risk of sepsis!
- Implants are only permitted for single use, not for reuse.

WARNING

Use of damaged or defective instruments/implants

- Risk of injury due to premature implant failure!
- Instruments/Implants with identifiable damage may not be used.
- Avoid notches, scratches or bending the implant in order to preserve its stability.

WARNING

Use of implants contrary to intended purpose

- Risk of injury due to implant failure!
- Implants must only be used in accordance with intended purpose.

WARNING

Combination of implant components of different sizes

- Damage to implant components!
- Combine only components of the same size.

WARNING

Use of instruments with electrical energy

- Risk of injury due to implant failure!
- Do not damage the surfaces of the implants under any circumstances.

Note: The following warning applies only to the foot plate systems.

WARNING

Bending to adjust

- Risk of injury due to implant failure!
- Do not overbend the implant.
- Do not bend implants repeatedly, as it weakens the implant material.
- Do not bend above the threaded holes.
- Using unsuitable instruments can lead to implant failure.

Implants must be inserted with the appropriate Merete instruments. The use of Merete instruments for any other purpose is not permitted.

The implants are to be used in a standardised operating theatre under appropriate hygiene and lighting conditions. The product may only be used by qualified surgeons who have substantial knowledge and experience in the field of orthopaedics and have demonstrated this proficiency accordingly.

Do not use any components that have been damaged during handling under any circumstances. When inserting and repositioning implants, the surgeon must ensure that implant surfaces have not been scratched or dented in any way. Even a tiny scratch can significantly reduce the life of an implant. The appropriate Merete instruments must be used for both implantation and explantation. The use of Merete instruments for purposes other than those intended is not permitted. Implants may not be reused under any circumstances.

Completed bone growth is a prerequisite for implant use.

Implantation in children:

- Only use the implants for children in exceptional cases. The decision is up to the user.
- Adjust the size and type of the implant as best as possible to the age, height, weight and bone development of the child.
- Due to growth, a subsequent operation for removal and/or revision is more likely in children.
- Perform follow-up checks at shorter intervals.

3.5.2. Sterile Implants

WARNING

Risk of infection due to non-sterile implants!

- Do not use implants whose packaging is damaged.
- Do not use implants beyond their expiry date.

WARNING

Use of contaminated instruments/implants

- Risk of sepsis!
- Use only instruments/implants without visible contamination.
- Handle instruments/implants only with sterile surgical gloves.

WARNING

Use of damaged or defective instruments/implants

- Risk of injury due to premature implant failure!
- Instruments/Implants with identifiable damage may not be used.
- Avoid notches, scratches or bending of the implant in order to preserve its stability.

Note: The following warning applies only to the foot plate systems.

WARNING

Implantation of trial implants

- Risk of injury due to breakage of trial implant!
- Only use trial implants in order to select a suitable permanent implant.
- Trial implants are not suitable for permanent implantation.

Implants which are delivered sterile are clearly labelled "STERILE". Sterile implants have been sterilised with 25 - 42 kGy (2.5 - 4.2 Mrad) gamma rays. All sterile products should be stored unopened in their original packaging in a cool – but frost free – dark, dry place until they are to be used. Before using the product, check the sterilisation expiration date on the product label, and check the protective packaging for damage. The red sterilisation dot on the package s/implants which are delivered sterile are clearly labelled "STERILE". Sterile implants have been sterilised with 25 - 42 kGy (2.5 - 4.2 Mrad) gamma rays. All sterile products should be stored unopened in their original packaging in a cool – but frost free – dark, dry place until they are to be used. Before using the product, check the sterilisation expiration date on the product label, and check the protective packaging for damage. The red sterilisation dot on the package serves as an indicator that the product is sterile.

Do not use products in damaged packaging. Observe the rules of asepsis when removing products from their protective packaging.

Generally speaking, implants delivered in trays or cartridges are not sterile.

3.5.3. Resterilisation

WARNING

Resterilization of implants

- Risk of injury due to premature implant failure caused by adverse material changes!
- Implants delivered sterile by Merete GmbH must not be resterilized and/or repacked.
- Products whose expiry date has passed may be returned to Merete GmbH.

Merete GmbH wishes to provide patients with cutting-edge technology products at all times. As such, we do not offer product resterilisation.

3.5.4 Storage and Handling

WARNING

Risk of infection due to non-sterile implants!

- Do not use implants whose packaging is damaged.
- Do not use implants beyond their expiry date.

WARNING

Use of contaminated implants

- Risk of sepsis!
- Use only instruments/implants without visible contamination.
- Only handle instruments/implants while wearing sterile surgical gloves.

Implants must be handled with the greatest of care. Even the slightest surface defects increase wear and tear, which could lead to complications. Therefore, do not mark the prosthesis surfaces in any way and avoid touching them with metallic or hard objects (e.g. instruments). Any mechanical manipulation or alteration of implants negates their approval for use and is not permitted. Do not implant any components of the prosthesis that have been handled incorrectly or show signs of damage (e.g. scratches).

3.6. Handling of Instruments and non-sterile Implants

3.6.1. Cleaning, disinfection and sterilisation instructions for reusable surgical instruments and non-sterile implants

Procedure

- Cleaning
- Disinfection
- Sterilisation with hot steam (DIN EN ISO 17665-1)

3.6.2. Warnings

Instruments/implants which are supplied non-sterile are clearly labelled "NON-STERILE". Non-sterile instruments/implants must be cleaned, disinfected (if applicable) and sterilised prior to use. Only qualified personnel may prepare instruments/implants for use. Use only approved cleaners and disinfectants (RKI, DGHM/VHA, FDA) (pH < 12 for metal instruments). Plastic instruments and instruments with plastic components must not be sterilised using dry heat under any circumstances. Instruments/implants supplied sterile are intended for single use only. Reuse or resterilisation is not permitted.

NOTE**Sterilization of instruments/implants supplied non-sterile**

If Merete products are sterilized by the user, this must be noted in the surgical report. All relevant labels and user instructions must be retained.

→ Observe current RKI guidelines.

→ Observe the standard preparation instructions provided.

3.6.3. Limitations on reprocessing

Other manufacturing information and recommendations must be followed.

3.6.4 Location of Use

It is advisable to process instruments for reuse as soon as possible after using them. Remove surface soiling with a disposable/paper towel. Immediately after use, the instrument may be laid in a disinfectant solution or in hot water 80°C (176°F), in order to reduce the risk of infection for the user and to facilitate cleaning.

3.6.5 Preparation for Cleaning/Pre-Cleaning

Dismantle all instruments/implants as far as possible. Take care not to lose any small components or screws. Pre-cleaning of instruments/implants: completely immerse the instrument in an enzymatic or alkaline cleaning solution (pH ≤ 12) and soak for 10 min. Clean the instrument/implant with a soft plastic brush. Then rinse the instrument/implant with running water for at least 1 minute.

Instruments/implants which cannot be disassembled, such as trial plates, forceps and side cutters, must be thoroughly cleaned by hand, with due attention to internal spaces (or cavities).

3.6.6 Automatic Cleaning and Disinfection

Mechanised cleaning is preferable to manual cleaning, if this is an option. Manual pre-cleaning improves automatic cleaning performance. The machine should offer a suitable thermal disinfection programme: A0 value > 3000, or at least 10 minutes at 93°C in older machines. Alternatively, if using a chemical disinfection method, bear in mind the risk of residue being left on the instruments. When choosing a cleanser, make sure that it is compatible with instrument materials. Follow manufacturer instructions when loading cleaning machines. Place instruments in such a way so as to allow complete, thorough rinsing of all ducts and cavities. Use deionised water for the final rinse. Be sure to allow sufficient drying time. Immediately after the programme has finished, remove instruments from the machine and, if necessary, dry them with a soft, absorbent, lint-free cloth.

3.6.7 Manual Cleaning

Pre-clean instruments manually to remove heavy surface contamination. Never use abrasive cleaning agents or metal brushes. Lay the instrument into the cleaning solution, following manufacturer's instructions as regards concentration, soaking time, and compatibility with the instrument materials. It must be ensured that the instrument is fully covered in cleaning solution and that all cavities, lumens and openings can be vented. Clean lumens and drill holes using appropriate brushes. After cleansing, rinse using deionised water, and dry thoroughly. For subsequent ultrasonic cleaning, make sure that the ultrasonic bath is preheated according to device manufacturers or cleaning agent manufacturer's instructions. When loading the bath, make sure that the cleaning solution completely covers the instruments and that all cavities, lumens and openings are fully vented. Cleaning of the instruments proceeds at 35-40 kHz for 5 minutes. After UV cleansing has finished, rinse instruments thoroughly with clear, running water, making sure to flush out cavities, lumens and openings wherever applicable.

3.6.8 Instrument Maintenance

Allow instruments to cool down to room temperature. Oil moving parts lightly with sterilisable, steam-penetrable surgical lubricating oil.

3.6.9 Checking Functionality

Use of contaminated instruments

- Risk of sepsis!
- Use only instruments without visible contamination.
- Handle instruments only with sterile surgical gloves.

Use of damaged or defective instruments/implants

- Risk of injury due to premature implant failure!
- Instruments/implants with identifiable damage may not be used.

NOTE

Sterilisation of instruments/implants supplied non-sterile

If Merete products are sterilised by the user, this must be noted in the surgical report. All relevant labels and user instructions must be retained.

- Observe current RKI guidelines.
- Observe the standard preparation instructions provided.

After each cleaning/disinfection, the instruments must be inspected for cleanliness, functionality and damage (e.g., bent, broken, worn or missing parts). Damaged instruments must be discarded and may not be used again. The instrument trays must be checked for completeness on the basis of the attached tray inserts. Training in the correct handling of implants and instruments by an authorised Merete representative is essential.

3.6.10 Cleaning and Disinfecting Empty Trays

Clean and disinfect empty trays using the same procedure, and under the same conditions, as for instruments. Be sure that the tray is completely dry before inserting instruments.

3.6.11 Packing

Prior to steam sterilisation, cleaned, disinfected instruments should be packed in suitable containers or sterilisation packs (DIN EN ISO 11607-1).

3.6.12 Sterilisation

Instrument/implant sterilisation is to be performed using an DIN EN ISO 17665 validated procedure. We recommend a fractionated vacuum method when using a steam steriliser. The maximum allowable sterilisation temperature is 134°C (273°F) plus tolerances as per DIN EN ISO 17665-1. Sterilisation time (exposure time at the sterilisation temperature) should be at 5 min at 134°C (273°F). The recommended pressure level is 3 bar.

3.6.13 Storage

After sterilisation, the instruments must be kept in their sterilisation packaging and stored in a dry, dust-free place.

3.6.14 Processing Instructions in accordance with DIN EN ISO 17664

The processor is responsible for ensuring that the processing procedure actually employed (using the materials, equipment, and personnel available in the processing facilities) achieves the desired results. Normally, this means that the procedure must be validated and subject to routine monitoring. Likewise, any deviation from the instructions provided should be carefully evaluated by the processor to determine its effectiveness and possible negative consequences. As regards the processing procedure, we would also like to point out the vital importance of following the Robert Koch Institute guidelines (published in the 2012, 55 issue of the German Federal Health Gazette) as well as any relevant national regulations.

3.6.15 Reusable Instruments

Handle instruments with care in order to ensure their long-term functionality. Unless Merete GmbH has specified otherwise, non-sterile implants/instruments may be reused indefinitely, provided that they are fully functional (see all points in 3.6.). Non-sterile implants/instruments with a defined limited period of usage must be labelled as follows: e.g. Q1/2019. On the expiry of the defined quarter, the implant/instrument must be inspected and approved by Merete GmbH to verify its continued functionality.

Operators are to use implants/instruments in accordance with RKI guidelines and MPBetreibV regulations.

4. Explantation

Explanation

- Risk due to bone loss!
- Increasing service life results in osseointegration; the surgeon/physician is responsible for decisions regarding implant removal.

5. Information Updates

Subject to technical change without notice.

Further information available from Merete GmbH or your authorised Merete GmbH representative.

New clinical experience may result in changes to certain procedures. Merete GmbH recommends regular attendance at refresher courses.

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1. Symbols

	CE mark		Do not reuse
	Manufacturer		Consult instructions for use
	Date of manufacture		Caution, consult accompanying documents
	Use-by date		Quantity
	Batch code		Hexalobe connection without cannulation
	Catalogue number		Hexalobe connection with cannulation
	Sterilized using ethylene oxide		Caution: Federal law (USA) restricts this device to sale by, or on order of a physician
	Sterilized using irradiation		
	Do not resterilize		
	Non-sterile		
	Do not use if package is damaged		
	Keep dry		
	Temperature limit		

Federal Law restricts this device to sale by or on the order of a physician.

2. General Information

Prior to using Merete products, surgeons and assisting staff are to study the safety information and product specific guidelines listed in these instructions for use as well as the surgical technique in detail. The relevant documentation is available from Merete on request. Surgeons must also be aware of any remaining risks associated with the products he or she intends to use and must inform patients of such risks in advance. Patients also need to be informed of any measures they can take to help prevent potential consequences of these risks. Implant operations must only be performed by surgeons who are not only qualified to carry out such operations, but also have extensive proven knowledge of, and experience in, this field. The surgeon bears all responsibility for negative consequences or complications arising from misdiagnosis, improper surgical technique, incorrect implant selection or handling, or failure to observe the safety instructions provided in this product information sheet. Neither the manufacturers nor authorized Merete product representatives may be held liable in such cases. Before operating, study the techniques outlined in the manual carefully.

Operation planning should be done on the basis of in-depth evaluations of patient X-rays, which provide the information necessary for determining the appropriate prosthesis type, size and possible combinations. The size of the osteosynthesis implant is selected based on the given anatomical conditions.

The implants are to be used in a sterile operating room.

2.1. Device description

All plates of the Locking Bone Plate Systems are made of Ti-6Al-4V ELI and have at least 2 K-wire holes for temporary bone fixation with K-wires.

Merete® MetaFix™ Small Fragment Locking Bone Plate System:

- MetaFix™ I Plates are T-shaped plates to be used with 4 MetaFix™ LS Locking Screws.

Merete® BLP™ Small Fragment Locking Bone Plate System:

- MetaFix™ BLP Plates are U-shaped plates to be used with 4 MetaFix™ LS Locking Screws.

Merete® Locking Bone Plate System I:

- Style 1 plates (MetaFix™ Ludloff) are U-shaped plates to be used with 2 MetaFix™ LS Locking Screws and 2 non-locking Merete® Cannulated PCS screws.
- Style 2 plates (MetaFix™ MTP) are H-shaped plates to be used with 6 or 7 MetaFix™ LS Locking Screws and either 1 non-locking Merete® Cannulated PCS Screw or 1 non-locking Merete® CS Cortical Screw.

Merete® Locking Bone Plate System II:

- Style 1 plates (MetaCun™) are straight plates to be used with 2 MetaFix™ LS Locking screws.
- Style 1 plates (MetaCun™ II) are straight plates to be used with 2 MetaFix™ LS Locking screws and 1 non-locking Merete® CS Cortical Screw.
- Style 2 plates (DuoMetaCun™) are U-shaped plates to be used with 4 MetaFix™ LS Locking screws.
- Style 2 plates (DuoMetaCun™ II) are U-shaped plates to be used with 4 MetaFix™ LS Locking screws and 2 non-locking Merete® CS Cortical Screws.
- Style 6 plates (MetaFix™ OpenWedge) are H-shaped plates to be used with 4 MetaFix™ LS Locking Screws.

Merete® Locking Bone Plate System III:

- Style 10 plates (MetaFix™ BG10) are L-shaped plates to be used with 4 MetaFix™ LS Locking Screws.
- Style 11 plates (MetaFix™ PlantarMAX™) are U-shaped plates to be used with 4 MetaFix™ LS Locking Screws and 1 non-locking Merete® Cannulated PCS Screw.
- Style 13 plates (MetaStep™ Calcaneal Step Plate) are stepped H-plates to be used with 6 Meta-Fix™ LS Locking Screws.
- Style 14 plates (SCARFix™) are L-shaped plates to be used with 2 MetaFix™ LS Locking Screws and 2 or 3 Merete® Cannulated PCS Screws.

MetaFix™ LS Locking Screws:

The screws are fully threaded and self-tapping with a threaded head to lock into Merete Locking plates. Locking screws/plates incorporate a screw-to-plate locking feature which creates a locked, fixed angle construction to hold fracture or osteotomy reduction. The screws are made of titanium alloy Ti-6Al-4V ELI and are available in the diameters 3.0 mm, 3.5 mm and 3.8 mm. 3.0 mm diameter is available in 12 mm to 32 mm in 2 mm increments; 3.5 mm diameters is available in 12 mm to 48mm in 2 mm increments excluding 42 mm and 46 mm; 3.8 mm diameters is available in 12 mm to 48 mm in 2 mm increments.

Merete® Cannulated PCS Screws:

The Merete® Cannulated PCS Screws are cannulated, self-drilling and self-tapping screws for the compression of MetaFix™ Plates as well as for the fixation of small bone fractures. The screws are made of titanium alloy Ti-6Al-4V ELI and are available in diameter 3.0 mm in the lengths from 12 mm to 40 mm to 48 mm in 2 mm increments. Characteristics include a short-threaded shaft with a non-threaded head and 3 cutting flutes.

Merete® CS Cortical Screws:

The Merete® CS Screw is a self-tapping compression screw which can be

used with the MetaFix™ System as well as for the fixation of small bone fractures. The screws are made of titanium alloy Ti-6Al-4V ELI and are available in diameter 3.0 mm and 3.5 mm in the lengths from 12 mm to 40 mm in 2 mm increments. Characteristics include a fully threaded shaft with a non-threaded head and 3 cutting flutes.

Merete® Cannulated HCS Screws :

Merete® Headless Compression Screws are cannulated, self-drilling/self-tapping, dual pitched screws with a threaded head which can be countersunk into the bone. The screws are provided in the diameter 3.0 mm in the lengths from 12 mm to 32 mm and 35 mm and 40 mm in 2 mm increments, in diameter 4.3 mm in the lengths from 24 mm to 50 mm and 55 mm in 2 mm increments and in diameter 6.5 mm in the lengths from 40 mm to 120 mm in 5 mm increments in various lengths and are manufactured from titanium alloy Ti-6Al-4V ELI.

Mecron® Cannulated Screw System:

The Mecron® Cannulated Screw System consists of self-tapping and self-drilling screws which are manufactured from titanium alloy Ti-6Al-4V ELI. The Mecron® Cannulated Screw System is comprised of headed and headless bone screws and corresponding washers for headed screws. The screws are offered in diameters from 2.0 mm thru 4.0 mm (in 0.5 mm diameter increments) and overall lengths from 8 mm (for smaller diameters) thru 50 mm in 2.0 mm increments.

Merete® TwistCut™ Snap-Off Bone Screw:

The self-drilling and self-tapping TwistCut™ Snap-Off Screw is made of titanium alloy Ti-6Al-4V ELI. The shaft separates from the screw when the head meets cortical bone. The screw is 2.0 mm in diameter and 11 mm to 14 mm in length.

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All Merete Locking bone plates are used in combination with Merete screws for adult and pediatric patients. Indications for use include fixation of fresh fractures, revision procedures, joint fusion and reconstruction of small bones of the hand, feet, wrist, ankles, fingers and toes.

The MetaFix™ LS Locking Screws are solely intended to be used in combination with Merete Locking plates for adult and pediatric patients as indicated above.

The Merete Compression Screws, Merete® CS Screw and Merete® Cannulated PCS Screws, can be used in combination with Merete Locking bone plates or individually. The Merete® Cannulated HCS Screws can only be used individually.

They are indicated for fracture fixation and reconstruction of various bones, including, osteotomies in the foot (as Hallux Valgus treatment) or hand, arthrodesis in hand, foot or ankle surgery, fixation of bone fragments in long bones or small bone fractures. The size of these screws should be adapted to the specific indication.

The Mecron® Cannulated Screw System is intended for use in bone reconstruction, osteotomy, arthrodesis, joint fusion, ligament fixation, fracture repair and fracture fixation appropriate for the size of the device.

The Merete® TwistCut™ Snap-Off Bone Screw is indicated for fixation of bone fractures or for bone reconstruction. Examples include Fixation of Small Bone Fragments, Weil-Osteotomy, Mono-Cortical Fixation, Osteotomies and fracture fixation in the foot and hand.

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- Risk of injury due to premature implant failure!
- Instruments/implants with identifiable damage may not be used.
- Avoid notches, scratches or bending of the instrument/implant in order to preserve its stability.

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- Damage to/destruction of instrument/implant and injury to patient!
- Ensure correct implant/instrument handling.
Do not misuse!

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WARNING

Improper use of an implant/instrument

- Damage to/destruction of instrument/implant and injury to patient!
- Ensure correct implant/instrument handling. Do not misuse.
- Ensure to use the instruments belonging to the implant.

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WARNING

Use of contaminated instruments/ implants

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WARNING

Risk of infection due to non-sterile implants

- Do not use implants whose packaging is damaged.
- Do not use implants beyond their expiration date.

WARNING

Reste­ri­za­tion of im­plants

- Risk of injury due to premature implant failure caused by adverse material changes!
- Implants delivered sterile by Merete Technologies, Inc. must not be resterilized and/or repacked.
- Products whose expiration date has passed may be returned to Merete Technologies, Inc.

NOTE

Sterilization of instruments/implants supplied non-sterile

- If Merete products are sterilized by the user, this must be noted in the surgical report. All relevant labels and user instructions must be retained.
- Observe the standard preparation instructions provided.

NOTE

Non-sterile screws, whose color cannot be clearly identified must not be used.

- The Merete Locking Screws are intended to be used in combination with Merete Locking Plates only.
- Plates and Screws are for single use only. An internal fixation device must never be reused.
- Non-sterile devices must be properly cleaned and sterilized before use.

2.5. Precautions

- Merete implants are designed for use by surgeons experienced in the appropriate specialized procedures. It is the responsibility of the surgeon to become familiar with the proper surgical techniques.
- Prior to use, examine device and check for proper functioning.
- The appropriate Merete instruments must be used for both implantation and explantation. The use of Merete instruments for purposes other than those intended is not permitted.
- Additional precautions include those applicable to any surgical procedure. In general, careful attention must be paid to asepsis and avoidance of anatomical hazards.
- Completed bone growth is a prerequisite for implant use. Only use the implants for adolescents in exceptional cases. The decision is up to the user. Adjust the size and type of the implant as best as possible to the age, height, weight and bone development of the adolescent. Due to growth, a subsequent operation for removal is more likely in adolescents. Perform follow-up checks at shorter intervals.
- For use in pediatric patients this product should not be placed across the growth plate.
- Explantation: Risk due to bone loss. Implant osseointegration may occur over time. The surgeon/ physician is responsible for decisions regarding implant removal.
- For the safe explantation of ingrown screws non-cannulated screw drivers (Ref. A114332 or Ref. A114333) with fixed handles are required. DO NOT use the cannulated screwdriver for explantation.

2.6. Factors Interfering with success

- Obesity and pre-obesity
- Local bone tumours
- Osteoporosis or osteomalacia
- Systemic or metabolic disorders
- Smoking, alcohol or drug abuse

- Physical activities associated with strong shocks which could result in the implant being exposed to impacts and/or excessive loads (e.g. hard physical labour, certain types of sport)
- Patients with pathological mental or neurological conditions, or patients who are not capable of following the necessary post-operative care instructions

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- Loosening and dislocation of the implant due to unfavourable anchoring/suboptimal positioning in the bone
- Implant failure due to excessive stress
- Early and late onset infections
- Wound haematomas and delayed wound healing
- Vascular disorders, including vein thrombosis, pulmonary embolism and heart failure
- Heterotopic ossification
- Nerve damage as a result of surgical trauma
- Difficult removal of the implant with long period of implantation (> 6 months) due to grown in bone
- Pain (chronic regional pain syndrome, hardware-related)
- Pseudarthrosis / Insufficient healing of bone fracture sites
- Delayed fracture healing
- Misaligned fracture healing
- Improper fixation
- Implant failure, implant breakage and/or metal wear
- Loss of correction
- Soft tissue irritation
- Pressure

- Scar infection and pain
- Incision dehiscence
- Neuritis
- Under- or overcorrection
- Uncontrolled fracture of the lateral corticalis
- Shortening of the first metatarsals
- High recurrence rate of hallux valgus in growing patients
- Transfer metatarsalgia
- Necrosis (avascular and soft tissue)
- Extremity Amputation
- Ulcer
- Skin Numbness and inflammation

Forefoot implants

- Impaired function of the metatarsophalangeal I joint
- Osteonecrosis of the first metatarsal head
- Arthrosis of the metatarsophalangeal I joint
- Lesion of the extensor or flexor tendon resulting in functional impairment
- Painful osteophyema of the metatarsophalangeal joint
- Contracture of the interphalangeal joint

Midfoot implants

- Tendinitis of the M. tibialis anterior
- Tendinitis of the M. extensor hallucis longus
- Capsulitis of the metatarsophalangeal II joint
- Sesamoidosis
- Deformity of the metatarsal I

- Residual hallux interphalangeus
- Rupture of the tibialis anterior tendon

Hindfoot implants

- Calcaneodynia
- Peroneal impingement syndrome
- Traumatic arthritis

2.8. MRI Safety Information

The Merete[®] Locking Bone Plate Systems, MetaFix[™] LS Locking Screws, Merete[®] Cannulated PCS, Merete[®] CS Cortical Screws, Merete[®] Cannulated HCS, TwistCut[™] Snap-Off and Mecron[®] Cannulated Screws have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the plates and screws in the MR environment is unknown. Scanning a patient who has these plates/screws may result in patient injury.

CAUTION

Examination of patient using MRI

- Risk of injury due to alternating magnetic fields!
- Merete Technologies, Inc. does not authorize the use of MRI examinations in conjunction with the components described in these user instructions.
- Always perform an individual risk-benefit analysis.
- Check whether other imaging procedures can be used to achieve the desired diagnostic goal.

3. Further information for use

Cleaning, disinfection, and sterilization instructions for reusable surgical instruments and non-sterile implants.

WARNING

Use of instruments contrary to intended purpose

- Risk of injury due to instrument failure!
- Instruments must only be used in accordance with intended purpose.

Process

- Cleaning
- Disinfection
- Sterilization with hot steam (DIN EN ISO 17665-1)

3.1. Warnings

WARNING

Use of damaged or defective instruments

- Risk of injury due to premature instrument failure!
- Instruments with identifiable damage may not be used.
- Avoid notches, scratches or bending of the instrument in order to preserve its stability.

WARNING

Use of instruments contrary to intended purpose

- Damage to/destruction of instrument and injury to patient!
- Ensure correct instrument handling. Do not misuse!

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- Risk of sepsis!
- Use only instruments/implants without identifiable soiling.
- Handle instruments/implants only with sterile surgical gloves.

WARNING

Risk of infection due to non-sterile instruments

- Do not use instruments whose packaging is damaged.
- Do not use instruments beyond their expiration date.

Instruments/implants that are supplied non-sterile are clearly marked with “NON STERILE” and must be cleaned, disinfected and sterilized prior to use. The instruments/implants may only be processed by qualified personnel. Only approved cleaning agents and disinfectants (refer to FDA Guidance) are to be used ($\text{pH} \leq 12$ pH for instruments/implants made of metal). Instruments made from synthetic materials or ones that contain components made from synthetic materials are not to be sterilized using dry heat.

3.2. Restrictions regarding reprocessing

The presented cleaning processes have been validated. Other methods of cleaning may be suitable but need to be validated by the user of the device. Differing instructions and recommendations of the manufacturer are to be observed.

NOTE

Non-sterile screws, whose color cannot be clearly identified must not be used.

3.3. Point-of-use processing

It is advisable to prepare instruments/implants for reprocessing as soon as possible after having used them. Macroscopic surface contamination can be difficult to remove by automated cleaning procedures. Prior to cleaning remove macroscopic contamination with a disposable towel/paper towel. Keep instruments/implants moist after use to prevent contamination from drying. Instruments/implants may be placed in a disinfectant solution or hot water 176°F (80°C) immediately after use in order to facilitate cleaning and to reduce risk of infection.

3.4. Preparation for cleaning

Dismantle all instruments as completely as possible. Be sure to ensure to keep all small components and screws. Pre-cleaning of instruments/implants:

- Completely immerse the instruments/implants in an enzymatic or alkaline cleaning solution ($\text{pH} \leq 12$) and soak for 10 min
- Clean the instruments/implants with a surgical scrub brush
- Then rinse the instrument/implant for at least 1 min with deionized water

3.5. Automatic cleaning

Automatic cleaning is preferable to manual cleaning, if this is an option. The machine should offer a suitable thermal disinfection program. Minimum cycle steps and parameters are:

- Rinse 1 min with cold water $< 109.4^\circ\text{F}$ ($< 43^\circ\text{C}$)
- Cleaning 5 min with cleaning agent (131°F (55°C) or follow the manufacturer's instructions)
- 1 min neutralization with warm water
- 1 min rinse (note: final rinsing is to be carried out with deionized water)

When choosing a cleanser, make sure that it is compatible with instrument/implant materials. Follow manufacturer's instructions when loading cleaning machines. Place instruments/implants in such a way so as to allow complete, thorough rinsing of all ducts and cavities. Use deionized water for the final rinse. Immediately after the cycle has completed, remove instruments/implants from the machine and, if necessary, dry them with a soft, absorbent, lint-free cloth. Ensure to allow sufficient drying time.

3.6. Automatic disinfection

Choose a program for an A0 value > 3000 , or at least 10 minutes at 200°F (93°C) in older machines. Alternatively, if using a chemical disinfection method, bear in mind the risk of residue being left on the instruments/implants.

3.7. Manual cleaning

Begin by removing major surface contamination from the instruments/implants using a soft nylon brush or a soft, lint-free cloth, along with either clear running water or a cleaning solution. Never use abrasive cleaning agents or metal brushes.

Place the instruments/implants in the cleaning solution, following manufacturer's instructions with regard to concentration, soaking time, and compatibility with the instrument/implant materials. Ensure that the instruments/implants are completely submerged in the cleaning solution. Be sure to vent all cavities, lumens and openings. Clean lumens and drill holes using appropriate brushes. After cleansing, rinse using deionized water, and dry thoroughly.

Subsequent ultrasonic cleaning: Be sure that the ultrasonic bath is pre-heated according to device manufacturer or cleaning-agent manufacturer's instructions. When loading the bath, make sure that the cleaning solution completely covers the instruments/implants, and that all cavities, lumens and openings are fully vented. Clean instruments/implants at 35-40 kHz for five minutes. After ultrasonic cleansing has finished, rinse instruments/implants thoroughly with deionized water, making sure to flush out cavities, lumens and openings wherever applicable.

3.8. Instruments maintenance

Allow instruments to cool down to room temperature. Lubricate moving parts lightly with sterilizable, steam-penetrable surgical lubricating oil.

3.9. Checking functionality

After each cleaning/disinfection, inspect the instruments for cleanliness, functionality, and damage (e.g., bent, broken, worn or missing parts). Never use damaged instruments or implants.

3.10. Cleaning and disinfecting empty trays

Clean and disinfect empty trays using the same procedure and under the same conditions as for instruments/implants. Be sure that the tray is completely dry prior to inserting instruments.

3.11. Packing

Prior to steam sterilization, cleaned, disinfected instruments/implants should be inserted into suitable containers or sterilization packages (DIN EN ISO 11607-1)

3.12. Sterilization

 <p>WARNING</p>	<p>Reste­ri­liza­tion of single use instru­ments</p> <ul style="list-style-type: none">– Risk of injury due to premature instrument failure caused by adverse material changes!→ Single use instruments delivered sterile by Merete Technologies, Inc. must not be resterilized and/or repackaged.→ Products whose expiration date has passed may be returned to Merete Technologies, Inc.
<p>NOTE</p>	<p>Steri­liza­tion of instru­ments/implants sup­plied non-sterile</p> <p>If Merete products are sterilized by the user, this must be noted in the surgical report. All relevant labels and user instructions must be retained.</p> <ul style="list-style-type: none">→ Observe the standard preparation instructions provided.

Implants which are delivered sterile are clearly labelled "STERILE". Sterile implants have been sterilised with 25 - 42 kGy (2.5 - 4.2 Mrad) gamma rays. All sterile products should be stored unopened in their original packaging in a cool – but frost free – dark, dry place until they are to be used. Before using the product, check the sterilisation expiration date on the product label, and check the protective packaging for damage. The red sterilisation dot on the package serves as an indicator that the product is sterile. Do not use products in damaged packaging. Observe the rules of asepsis when removing products from their protective packaging. Generally speaking, implants delivered in trays or cartridges are not sterile.

Sterile implants are delivered in sterile packing, sterilized by Gamma Radiation:

Method	Gamma Radiation
Radiation Dose	2.5 – 4.2 Mrad
Sterility Assurance Level	10^{-6}
Sterility Validation Method	ISO 11137-2 / $VD_{max}25$
Packaging	Double blister packaging: PETG, Tyvek
Pyrogenicity	Not labeled pyrogen free

Non-sterile implants and instruments are delivered non-sterile. Following sterilization method is recommended:

Method	Steam Sterilization
Cycle	Pre-Vacuum
Temperature	270°F [132°C]
Expose Time	4 minutes
Drying Time	20-30 minutes
Sterility Assurance Level (SAL)	10 ⁻⁶
Sterility Validation Method	Overkill method per ISO 17665-1:2006
Wrapping	FDA cleared wrap

The trays containing implants and instruments are only provided for packaging, not for sterilization.

3.13. Storage

WARNING

Use of contaminated instruments

- Risk of injury due to premature instrument failure!
- Use only instruments without visible contamination.
- Handle instruments only with sterile surgical gloves.

WARNING

Risk of infection due to non-sterile instruments

- Do not use instruments whose packaging is damaged.
- Do not use instruments whose expiration date has passed.

WARNING

Use of damaged or defective instruments

- Risk of injury due to premature instrument failure!
- Instruments with identifiable damage may not be used.
- Avoid notches, scratches or bending of the instrument in order to preserve its stability.

After sterilization, the instruments/implants must be kept in their sterilization packaging and stored in a dry, dust-free place.

3.14. Preparation instructions in accordance with DIN EN ISO 17664

The preparer is responsible for ensuring that the preparation procedure actually employed (using the materials, equipment, and personnel available in the preparation facilities) achieves the desired results. Normally, this means that the procedure must be validated and subject to routine monitoring. Likewise, any deviation from the instructions provided should be carefully evaluated by the preparer to determine its effectiveness and possible negative consequences.

3.15. Disposal

Disposal of implants according to standard hospital procedures.

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Merete Technologies Inc.