Important Product Information OsteoBridge™ IDSF (Package Insert OsteoBridge™ IDSF - Intramedullary Diaphyseal Segmental Defect Fixation System)**Caution:**

Federal law restricts this device to sale by or on the order of a physician.

1. Explanations to used symbols

	Manufacturer
	Use-by date
	Batch code
	Catalogue number
	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Do not reuse
	Consult instructions for use
	Caution
	Quantity
	Caution: Federal law (USA) restricts this device to sale by or on order of a physician
	Keep away from sunlight
	MR Conditional
	Single sterile barrier system with protective packaging outside
	Unique device identifier

2. Description

Warning: Use of damaged or defective implants. Risk of injury due to premature implant failure! Implants with identifiable damage may not be used. Avoid notches, scratches or bending of the implant in order to preserve its stability.

Warning: Use of damaged or defective instruments. Risk of injury due to premature implant failure! Instruments with identifiable damage may not be used.

Warning: Use of implants which have been previously used. Risk of injury due to premature implant failure! Risk of sepsis! Implants are only approved for single use, not reuse.

The OsteoBridge™ IDSF - Intramedullary Diaphyseal Segmental Defect Fixation System is intended to be used in the management of segmental diaphyseal bone loss secondary to radical bone loss and/or resection due to tumors in either humerus, tibia or femur in oncology patients. The device is divided into three categories/product size ranges that are suitable for each bony application (see Table 1, Table 2 and Table 3 below).

The system is designed for non-cemented and cemented fixation in bone. The modular system includes a spacer made of two cylindrical half-shells that clamp over the ends of two intramedullary nails with eight screws. Non-cemented implants are additionally fixed with two interlocking screws per nail. Spacers and intramedullary nails of different sizes can be selected to accommodate varying intramedullary canal diameters and different bone defect sizes. If the defect size is larger than 70 mm, two spacers are to be connected with a spacer connector. All components of the OsteoBridge™ IDSF Implant System are manufactured from Ti6Al4V ELI titanium alloy according to ASTM F136 and ISO 5832-3.

Table 1: Dimension for the OsteoBridge™ IDSF Humerus System

Humerus				
Spacer			Nail	
Outer DIA. [mm]	Clamping DIA. [mm]	Length [mm]	Nail DIA. [mm]	Length [mm]
20	10	40, 50, 60, 70	7	60, 70, 90
			8	60, 70, 90
			9	60, 70, 90
			10	60, 70, 90, 110, 130

Table 2: Dimensions for the OsteoBridge™ IDSF Tibia system

Tibia				
Spacer			Nail	
Outer DIA. [mm]	Clamping DIA. [mm]	Length [mm]	Nail DIA. [mm]	Length [mm]
25	14	40, 50, 60, 70	9	60, 70, 90
			10	60, 70, 90, 110, 130
			12	70, 90, 110, 130, 150
			14	90, 110, 130, 150, 200

Table 3: Dimensions for the OsteoBridge™ IDSF Femur system

Femur				
Spacer			Nail	
Outer DIA. [mm]	Clamping DIA. [mm]	Length [mm]	Nail DIA. [mm]	Length [mm]
34	16	40, 50, 60, 70	10	60, 70, 90, 110, 130
			12	70, 90, 110, 130, 150
			14	90, 110, 130, 150, 200
			16	110, 130, 150, 200
			18	90, 110
			20	90, 110

2.1. Intended Purpose/Use

Warning: Use of implants contrary to intended purpose. Risk of injury due to implant failure! Implant must only be used in accordance with intended purpose.

The OsteoBridge™ IDSF Implant System is intended to be used in the management of segmental diaphyseal bone loss secondary to radical bone loss and/or resection due to tumors in either humerus or tibia or femur in oncology patients. The modular system includes a spacer made of two cylindrical half-shells that clamp over the ends of two intramedullary nails with eight screws. Spacers and intramedullary nails of different sizes can be selected to accommodate varying intramedullary canal diameters and different bone defect sizes. If the defect size is larger than 70 mm, two spacers may be connected with a spacer connector. The nails can be used non-cemented and cemented. All components of the OsteoBridge™ IDSF System are manufactured from Ti6Al4V ELI titanium alloy according to ASTM F 136.

2.1.1. Indications**The system is indicated for**

- Long-term stabilization of major bone defects resulting from
 - Radical bone loss due to tumors and/or metastases
 - Bone resections following tumors and/or metastases
- For use only in the diaphyseal region of humerus, tibia and femur

2.1.2. Intended patient target groups

- Age: when the bone growth is completed
- Body weight: obesity or pre-obesity can interfere with the success of the implant
- Activity: observe aftercare, physical activities associated with strong shocks which could result in the implant being exposed to impacts and/or excessive loads (e.g., hard physical labour, certain types of sport)

2.2. Contraindications

- Profound osteoporosis
- Acute or chronic, local or systemic infections
- Any disease that can negatively influence the success of the operation
- Defective bone structures which would impede adequate anchoring of the implant
- Degenerative bone diseases
- Use in non-diaphyseal bone regions
- Adiposity

2.3. Warnings

Warning: Improper use of an implant/instrument. Damage to destruction of instrument/implant and injury to patient! Ensure correct handling of implant/instrument. Do not misuse. Ensure to use the instruments belonging to the implant.

Warning: Combination with products from other manufacturers. Risk of injury due to implant failure! Do not combine the implant components with products from other manufacturers.

Warning: Combination of implant components of different sizes. Damage to implant components! Combine only components of the same size.

Warning: Implantation of trial implants. Risk of injury due to failure of trial implant! Only use trial implants in order to select a suitable permanent implant. Trial implants are not suitable for permanent implantation.

Warning: Incorrect preparation of the medullary cavity. Risk of implant failure (lengthening/shortening of the extremity)! Drill into the medullary cavity as specified in the surgical technique and use a cement restrictor when using bone cement to ensure a proper cement mantle.

Warning: Improper use of bone cement. Risk of injury due to premature implant failure! Pay attention to the information given by the cement manufacturer. The cement mantle should be spread evenly 2mm on each side.

Warning: Foreign bodies (e.g. cement residues, tissue, bones) between implant components. Risk of injury due to implant failure! Thoroughly clean any foreign bodies from implant components.

Warning: Incorrect or improper assembly of implant. Risk of increased wear to implant! Risk of implant failure! Perform the implantation with care. Observe correct handling of implant components and instruments. A bicortical anchoring is essential for a stable implant when using interlocking screws. The interlocking screws must be completely tightened. Make sure there is at least 4 cm of cortical bone contact of the nail within the bone, in order to ensure stable anchorage. Prior to assembling the spacer, wash the clamping surface in order to completely remove all kinds of debris, including bone splinters, soft tissue parts, bone cement and others. It is essential to ensure that the nails of the implants are clamped into the spacer according to their respective minimal clamping length. The instructions for

the assembling of the spacer with the spacer screws are to be followed exactly. All eight screws must be used for spacer assembly. Use the guiding pins during assembly to ensure the spacer half-shells remain parallel to each other. When screwing in the guiding pins or spacer screws, be careful not to damage the spacer threads. In case of defective bone structure, the nail must be long enough to bridge over the damaged region.

Warning: Torque limiter not used or incorrectly used during implant assembly. Risk of increased wear to implant! Risk of implant failure! Use the torque limiter during assembly and follow the sequence instructed for tightening the screws to ensure even and sufficient clamping of the spacers. Tighten all eight screws three times in the given order until the correct torque has been reached. When reaching the correct torque a "pop" sound is triggered.

Warning: Use of instruments with electrical energy. Risk of injury due to implant failure! Do not damage the surfaces of the implants under any circumstances.

- This product is for single use only. An internal fixation device must never be reused. The reuse of explanted components is not permitted under any circumstances.
- During the implantation and repositioning, the operating surgeon should pay attention that the surfaces of implant are not damaged due to nicks and scratches. Even a slight scratch can considerably reduce the lifespan of an implant.

2.4. Cautions

- Physical activities or trauma associated with strong shocks which could result in the implant being exposed to impacts and/or excessive loads or sports activities that affect the affected limb or cause a component to fail.
- The devices are designed for use by surgeons experienced in the appropriate specialized procedures. It is the responsibility of the surgeon to read the surgical technique handbook as well as become familiar with the correct use of the instrumentation for this implant.
- Training in the correct handling of implants and instruments is only to be executed by an authorized Merete representative.
- Additional precautions include those applicable to any surgical procedure. In general, careful attention must be paid to asepsis and avoidance of anatomical hazards.
- The implants must be stored in the original packaging in a dry, clean place and at a temperature ranging from room temperature to cool but above frost-formation.
- Non-sterile devices must be properly cleaned and sterilized before use.
- The implantation must be performed with the corresponding Merete instruments. The use of Merete instruments for other than the intended purposes is not permitted.
- When implanting OsteoBridge IDSF the correlation tables are to be followed. The tables illustrate the possible combinations of the components with each other.
- Use trial components during the operation in order to check correct implants size and positioning.
- When using the spacer clamp, it is vitally important to ensure that the spacer half-shells are anchored in the correct position. This is the only way of ensuring that the spacer will be installed correctly.

Implantation in adolescents:

- Complete bone growth is a prerequisite for the application of the OsteoBridge™ IDSF System.
- Only use the OsteoBridge™ IDSF System for adolescents in exceptional cases. The decision is up to the user.
- Adjust the size and type of the implant as best as possible to the age, height, weight and bone development of the adolescents.
- Due to growth, a subsequent operation for removal and/or revision is more likely in adolescents.
- Perform follow-up checks at shorter intervals.

CAUTION: (Not recommended) Interlocking of nails with cementation leads to:

- Loosening or breakage of the locking screw
- Loosening or breakage of the nail
- Breakage of the cement mantle
- Heat generation
- Breakage of the drill (drill in the instrument tray is not designed for use in bone cement) can occur.

NOTE For spanning large defects not more than two spacers should be joined together with a connector.

Prior to using Merete products, surgeons and assisting staff are to study the safety information and product-specific guidelines listed in these instructions for use as well as the surgical technique (OsteoBridge™ IDSF) in detail. The respective documents are provided by Merete. Surgeons must also be aware of any remaining risks associated with the products he or she intends to use and must inform patients of such risks in advance. Only specifically qualified surgeons who have profound knowledge and expertise in the field of bone and joint surgery may perform the implantation. The surgeon bears all responsibility for adverse effects or complications arising from misdiagnosis, improper surgical technique, incorrect implant selection or handling, or failure to observe the safety instructions provided in the product information sheet. Neither the manufacturer nor authorized Merete product representatives may be held liable in such cases. Before operating, study the techniques outlined in the manual carefully. The patient has to be informed about alternative methods of treatment and other systems with adequate intended uses. All implants represent individual components of a single system. This necessitates those surgeons observe particular care in ensuring that all system components are mutually compatible and used for the applications intended. Such products may only be used in combination with other Merete products belonging to the same system. For further information, consult the catalogues and surgical techniques for the products in question. In case of doubt, contact your sales representative or Customer Service/Sales at Merete Technologies Inc. (MTI) prior to operating. Implantations must be performed using the appropriate Merete instruments. Usage of Merete instruments for purposes other than those originally intended is prohibited. Under no circumstances may implants be reused. Products which are delivered sterile may not be resterilized. Individual OsteoBridge™ IDSF implant components are attached to one another using secure clamp connections. This is only possible when the components are precisely coordinated in terms of size. Acceptable combinations are indicated in the surgical technique handbook and in the product brochure. There are no known interactions with any medications.

General Information on Implants

- The use of OsteoBridge™ IDSF implants should only be considered if, upon careful assessment, they appear to represent a better solution than other bone-preserving options.
- Care of a severely deformed, unhealthy joint or bone using an implant can reduce pain, increase joint load-bearing capacities and help improve patient mobility.
- Over time, there is always a possibility that even stable implants may loosen for any of several reasons which could make a revision necessary.
- The OsteoBridge™ IDSF Implant System may be inserted non-cemented or cemented (PMMA) according to the specifications supplied by the cement manufacturer. Do not additionally lock the nail after cementation.

Preoperative Planning

Surgery planning should be done on the basis of in-depth evaluations of patient X-rays which provide the information necessary for determining the appropriate implant type, size, and possible combinations. If desired, X-ray templates for preoperative planning are available from Merete. It should be examined preoperatively whether the existing implant sizes are suitable for the patient.

It is also important to carry out preoperative tests on the patient to rule out allergic reaction to the implant materials. Use trial components during the operation in order to check correct implants size and positioning. The components planned for implantation must be on hand in all available sizes. It is absolutely essential to check the implant label (type, item number, material and size) against the information on the package before insertion. Use the additional enclosed labels for documentation in surgery reports and for re-ordering. Check that all implant surfaces are free of soiling.

Important: Always check sterilization expiration dates! Do not use if expiration date has passed.

Application Environment

The implants and instruments are to be used in a standardized operating room. The products may only be used by qualified surgeons having a profound knowledge and experience in the field of orthopedics.

Revision Surgery

If the revision surgery involves replacement of all the implant components, the same instructions apply to the treatment of the implants as for a first implantation. A partial change of implants leads to the mixing of old, pre-damaged components and new implants. In such cases, it is important that the component mating surfaces be in excellent condition. It is always recommended to completely replace all components.

NOTE For explanation, please order the instrument tray GB90010 and GB90011 plus the extractor GB90203.

Storage and Handling of Implants

OsteoBridge™ IDSF implants must be handled with extreme care, since even minimal damage of the surfaces can lead to a failure of the implant. Therefore, even the subsequent marking of prostheses surfaces or touching them with metal or other hand-objects (e.g. instruments) must be avoided. The mechanical processing and modification of medical components leads to a loss of authorization and is not permitted. Implants that were handled in an improper way or exhibit damage (e.g. scratches) must not be implanted.

Patient Information

The patient must be informed by the doctor of the information provided in this document regarding factors which could interfere with the success of a surgery as well as possible complications which could result from such

an indication. The patient must be informed about which measures they can take to reduce the possible effects of these factors. All information given to the patient must be documented by the operating surgeon in writing.

All metal components are made of the following material:

- TiAl6V4 ELI alloy (ASTM F 136, DIN EN ISO 5832-3)

Further information about the chemical composition and the mechanical properties of the used materials can be requested from Merete.

2.5. Factors which could impair the Success of the Operation

- Allergies to the implant materials (particularly metals, but also bone cement)
- Inadequate bone density
- Systemic or metabolic disorders
- Smoking, alcohol or drug abuse
- A strong curvature of the medullary canal increases the risk of fractures during insertion of straight nails into bone
- When utilized at the humerus, no loads greater than 11 lb (5kg) should be lifted with the angled arm and
- Patients with mental or neurological illness or patients who are unable to follow necessary postoperative treatments.

2.6. Potential Adverse Effects

The potential adverse effects associated with this device are the same as with any metallic internal fixation device. These include but are not limited to the following:

- Loosening of the implant resulting from changed conditions of the load transfer or wearing and fracture of the cement mantle and/or tissue reaction of the implant
- Bending or fracture of the implant
- Early and late infection
- Undesired shortening, lengthening or malalignment of the relevant extremity as a result of suboptimal positioning of the implant
- Bone fractures, resulting from unilateral stress or weakened bone substance, reduction of bone density as a result of stress shielding or bone resorption as tissue reaction to abrasion articles
- Temporary or chronic neural damage resulting from pressure or hematoma
- Wound hematoma and delayed wound healing
- Vascular disease including venal thrombosis, pulmonary embolism and cardiac arrest
- Heterotopic ossification
- Nerve damage due to surgical trauma
- Metal sensitivity or allergic reaction to a foreign body
- Pain, discomfort, or abnormal sensation due to the presence of the device
- Delayed or non-union which may lead to breakage of the implant
- For non-cemented implantation lack of bony ongrowth onto the nails can lead to overloading of the locking screws

3. MRI Safety Information

MR Conditional
MR Safety Information/Indications for Use
Non-clinical testing has demonstrated that the Merete OsteoBridge™ IDSF (Intramedullary Diaphyseal Segmental Defect Fixation) System (consisting of two cemented or non-cemented intramedullary stems and one spacer or two spacers made of TiAl6V4 ELI [ISO 5832-3]) is MR conditional.
A patient with the entire assembled Merete OsteoBridge™ IDSF (Intramedullary Diaphyseal Segmental Defect Fixation) System can be safely scanned in an MR system meeting the following conditions:
- Static magnetic field of 1.5 Tesla and 3.0 Tesla
- Maximum spatial gradient field of 3,000 Gauss/cm (30 T/m).
- Maximum MR system reported whole-body-averaged specific absorption rate (SAR) at 1.5 Tesla or 3 Tesla of 1 W/kg for 15 minutes of scanning. Under the scan conditions defined above, the Merete OsteoBridge™ IDSF (Intramedullary Diaphyseal Segmental Defect Fixation) System is expected to produce a maximum temperature rise of less than 6 °C after 15 minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by the Merete OsteoBridge™ IDSF (Intramedullary Diaphyseal Segmental Defect Fixation) System extends at least 1 cm and up to approximately 7 cm from the device and exhibits substantial geometric distortion in the image when imaged with a gradient echo pulse sequence or a fast-spin echo pulse sequence and a 1.5 Tesla MRI system or a 3.0 Tesla MRI system.

4. Instrumentation

Please refer to instrumentation instructions for use (PT10027) for information regarding cleaning, disinfection and sterilization. All instruments are supplied non-sterile and must be sterilized prior to use according to the following recommended sterilization parameters.

Recommended sterilization parameters

Method	Temperature	Time	Drying Time
Steam Sterilization	270°F [132°C]	4 minutes	30 minutes

4.1. Process

- Cleaning
- Disinfection
- Sterilization with hot steam [DIN EN ISO 17665-1]

4.2. Warnings

Warning: Use of damaged or defective instruments. Risk of injury due to premature implant failure! Instruments with identifiable damage may not be used.

Warning: Use of contaminated instruments. Risk of Sepsis! Use only instruments without visible contamination.

Instruments that are supplied non-sterile are clearly marked "NON STERILE" and must be cleaned, disinfected and sterilized prior to use. The instruments may only be processed by qualified personnel. Only approved cleaning agents and disinfectants (FDA are to be used [pH ≤ 12 for metal instruments]). Instruments made from synthetic materials are not to be sterilized using dry heat.

4.3. Restrictions regarding Reprocessing

Warning: Resterilization of implants. Risk of injury due to premature implant failure caused by adverse material changes! Implants delivered sterile by Merete* GmbH must not be resterilized and/or repacked. Products or which expiry date has passed may be returned to Merete Technologies Inc. (MTI).

Warning: Use of damaged or defective instruments. Risk of injury due to premature implant failure! Instruments with identifiable damage may not be used.

The presented cleaning processes have been validated. Other methods of cleaning may be suitable but need to be validated by the user of the device. Instructions and recommendations of the manufacturer are to be observed.

4.4. Point-of-use-Processing

It is advisable to prepare instruments for reuse as soon as possible after having used them. Macroscopic surface contamination can be difficult to remove by automated cleaning procedures. Prior to cleaning, remove macroscopic contamination with disposable towel.

Keep instruments moist after use to prevent contamination from drying. Instruments may be placed in a disinfectant solution or hot water 176°F (80°C) immediately after use in order to facilitate cleaning and to reduce risk of infection.

4.5. Preparation for Cleaning

- Disassemble all instruments as far as possible. Ensure to keep all small components and screws. Pre-cleaning of instruments:
- Completely immerse the instruments in an enzymatic or alkaline cleaning solution (pH ≤ 12) and soak for 10 min
- Clean the instruments with a surgical scrub brush
- Then rinse the instrument for at least 1 min with deionized water

4.6. Automatic Cleaning

Automatic cleaning is preferable to manual cleaning, if this is an option. The machine should offer a suitable thermal disinfection program. Minimum cycle steps and parameters are:

- Rinse 1 min with cold water < 109.4°F (<43°C)
- Cleaning 5 min with cleaning agent [131°F (55°C) or follow the manufacturer's instructions]
- 1 min neutralization with warm water
- 1 min rinse (note: final rinsing is to be carried out with deionized water)

4.7. Automatic Disinfection

Choose a program for an A0 value > 3000, or at least 10 min at 200°F (93°C) in older machines. Alternatively, if using a chemical disinfection method, bear in mind the risk of residue being left on the instruments.

4.8. Manual Cleaning

Begin by removing major surface contamination from the instruments using a soft nylon brush or a soft, lint-free cloth, along with either clear running water or a cleaning solution. Never use abrasive cleaning agents or metal brushes. Place the instruments in the cleaning solution, following manufacturer's instructions regarding concentration, soaking time, and compatibility with all instrument materials. Ensure that the instruments are completely submerged in the cleaning solution. Be sure to vent all cavities, lumens and openings. Clean lumens and drill holes using appropriate brushes. After cleansing, rinse using deionized water, and dry thoroughly. For subsequent ultrasonic cleaning, be sure that the ultrasonic bath is preheated according to device manufacturer or cleaning-agent manufacturer instructions. When loading the bath, make sure that the cleaning solution completely covers the instruments, and that all cavities, lumens and openings are fully vented. Clean instruments at 35-40 kHz for 5 min. After ultrasonic cleansing has finished, rinse instruments thoroughly with deionized water, making sure to flush out cavities, lumens and openings wherever applicable.

4.9. Instruments Maintenance

Allow instruments to cool down to room temperature. Lubricate moving parts lightly with sterilizable, steam-penetrable surgical lubricating oil.

4.10. Checking Functionality

After each cleaning / disinfection, inspect the instruments for cleanliness, functionality and damage (e.g., bent, broken, worn or missing parts). Never use damaged instruments or implants. The completeness of instrument trays and cases should be assessed using the provided tray insert sheets.

4.11. Cleaning and Disinfecting Empty Trays

Clean and disinfect empty trays using the same procedure and under the same conditions as for instruments. Be sure that the tray is completely dry prior to inserting instruments.

4.12. Packing

Prior to steam sterilization, cleaned, disinfected instruments should be inserted into suitable containers or sterilization packages (DIN EN ISO 11607-1)

5. Sterilization

Warning: Risk of infection due to non-sterile implants! Do not use implants whose packaging is damaged. Do not use implants for whose expiry date has passed.

Warning: Use of contaminated implants/instruments. Risk of Sepsis! Use only implants/instruments without visible contamination.

Warning: Use of damaged or defective implants. Risk of injury due to premature implant failure! Implants with identifiable damage may not be used. Avoid notches, scratches or bending of the implant in order to preserve its stability.

Warning: Resterilization of implants. Risk of injury due to premature implant failure caused by adverse material changes! Implants delivered sterile by Merete* GmbH must not be resterilized and/or repacked. Products for which expiry date has passed may be returned to Merete Technologies Inc. (MTI)

NOTE Observe symbol on packaging: "Do not reuse".

OsteoBridge™ IDSF Nails, Spacers, Spacer Connectors and Interlocking Screws are delivered in sterile packaging, sterilized by Gamma Radiation:	
Method	Gamma Radiation
Radiation Dose	2.5 Mrad
Sterility Assurance Level	10 ⁻⁶
Sterility Validation Method	ISO 11137-2/VD _{max} 25
Packaging	Triple Peel PE Pouch (Nails) Double Peel PE Pouch (Other components)
Pyrogenicity	Not labeled pyrogen free

Instruments are delivered non-sterile. Following sterilization method is recommended	
Method	Steam Sterilization
Temperature	270°F [132°C]
Expose Time	4 minutes
Drying Time	20 - 30 minutes
Sterility Assurance Level (SAL)	10 ⁻⁶
Sterility Validation Method	Overkill method per ISO 17665-1:2006
Wrapping	FDA cleared wrap

The trays containing instruments are only provided for packaging, not to maintain sterility.

6. Storage

After sterilization, the instruments must be kept in their sterilization packaging and stored in a dry, dust-free environment.

The implants must be stored in the original packaging in a dry, clean place and at a temperature ranging from room temperature to cool but above frost-formation.

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