


Five advantages of the new BioBall[®] Instrument Tray



The top five advantages of the BioBall[®] Instrument Tray (HM30770)

- 1 **BioBall[®] AdapterSelector™** for tapers 12/14 and 14/16 – the first instrument for inspecting taper geometry intraoperatively.
- 2 **BioBall[®] Offset PositionAssistant** – for better visualisation of the offset adapter position ex situ.
- 3 New radiopaque **BioBall[®] Trial** components for visualisation via imaging techniques, inc. radiopaque BioBall[®] Trial Adapter 14/16 Offset.
- 4 **BioBall[®] Trial Adapters** are press-fit onto the stem taper.
- 5 Separating Instrument with adapter sleeve for **BioBall[®] Adapter** 12/14 S and 14/16 M.

Update
now

Why is there a new BioBall[®] Instrument Tray?

- All instruments and trial components are CE marked.
- Laser marking/printing for placement on the base of the tray.
- New clear and safe arrangement.
- Colour-coded holders for the test components.

Ref.	Description
HM30770	BioBall [®] Instrument Tray 12/14 and 14/16 non-steril