

## Discover the refined Merete BioBall® 12/14 Instrument Tray!

The Merete BioBall® Adapter System allows intraoperative adjustments without removing the well-fixed stem.

**THE POWER**

**TO ADAPT**

# Your **benefits** at a glance!

1

## **Practical Stackability**

The connecting feet on the tray base and lid enable secure, space-saving stacking of trays of the same size.



2

## **Enhanced Latching Mechanism**

The upgraded latches enable smooth, effortless handling, while rounded edges enhance safety and prevent damage to sterilization wraps, improving usability.



3

## **Durable Brackets**

The new nylon-coated aluminium brackets resist wear and tear, ensuring a robust and reliable tray design.



**THE POWER**

**TO**


**ADAPT**

# Order Information

## BioBall® Instrument Tray 12/14, Ref. HM30776

Included in Instrument Tray:

Ref.	Product Description	QTY
HI39006	BioBall® AdapterSelector® 12/14, non-sterile	1
HM10004	BioBall® Head Impactor, non-sterile	1
HM10005	BioBall® Universal Handle, non-sterile	1
HM10007	BioBall® Head Separating Wedge S, non-sterile	1
HM10008	BioBall® Head Separating Wedge M, non-sterile	1
HM10009	BioBall® Head Separating Wedge L, non-sterile	1
HM20001	BioBall® Adapter Extractor, non-sterile	1
HM20002	BioBall® Adapter Extractor Sleeve, for Adapter 12/14 S and MSPC M-L, non-sterile	1
HM39106	BioBall® Offset PositionAssistant 12/14, non-sterile	1
HM40028	BioBall® Trial Head Ø 28 mm, non-sterile	1
HM40132	BioBall® Trial Head Ø 32 mm, non-sterile	1
HM40036	BioBall® Trial Head Ø 36 mm, non-sterile	1
HM40040	BioBall® Trial Head Ø 40 mm, non-sterile	1
HM40121	BioBall® Trial Adapter 12/14, Standard, S (-3.0 mm), 12/14 (5°42'), non-sterile	1
HM40122	BioBall® Trial Adapter 12/14, Standard, M (0.0 mm), 12/14 (5°42'), non-sterile	1
HM40123	BioBall® Trial Adapter 12/14, Standard, L (+3.5 mm), 12/14 (5°42'), non-sterile	1
HM40124	BioBall® Trial Adapter 12/14, Standard, XL (+7.0 mm), 12/14 (5°42'), non-sterile	1
HM40125	BioBall® Trial Adapter 12/14, Standard, 2XL (+10.5 mm), 12/14 (5°42'), non-sterile	1
HM40126	BioBall® Trial Adapter 12/14, Standard, 3XL (+14.0 mm), 12/14 (5°42'), non-sterile	1
HM40127	BioBall® Trial Adapter 12/14, Standard, 4XL (+17.5 mm), 12/14 (5°42'), non-sterile	1
HM40128	BioBall® Trial Adapter 12/14, Standard, 5XL (+21.0 mm), 12/14 (5°42'), non-sterile	1
HM40222	BioBall® Trial Adapter 12/14, Offset, M (0.0 mm), 12/14 (5°42'), non-sterile	1
HM40223	BioBall® Trial Adapter 12/14, Offset, L (+3.5 mm), 12/14 (5°42'), non-sterile	1
HM40224	BioBall® Trial Adapter 12/14, Offset, XL (+7.0 mm), 12/14 (5°42'), non-sterile	1
HM40225	BioBall® Trial Adapter 12/14, Offset, 2XL (+10.5 mm), 12/14 (5°42'), non-sterile	1
HM40226	BioBall® Trial Adapter 12/14, Offset, 3XL (+14.0 mm), 12/14 (5°42'), non-sterile	1
HM40227	BioBall® Trial Adapter 12/14, Offset, 4XL (+17.5 mm), 12/14 (5°42'), non-sterile	1
HM40228	BioBall® Trial Adapter 12/14, Offset, 5XL (+21.0 mm), 12/14 (5°42'), non-sterile	1




**Order Information  
for 10 additional BioBall®  
Special Adapters**



**More Information  
on BioBall® Implants**

# Merete® Learning Center

Your central hub for extensive  
expertise and continued  
education in medical technology



[merete.de/en/learning-center](https://merete.de/en/learning-center)

## Merete GmbH

Alt-Lankwitz 102  
12247 Berlin  
Germany

Tel. +49 (0)30 77 99 80-0  
Fax +49 (0)30 76 68 03 61

[service@merete.de](mailto:service@merete.de)  
[merete.de/en](https://merete.de/en)



All content is protected by copyright and intellectual property rights owned by or licensed to Merete GmbH. Redistribution or duplication without written consent is prohibited. This document is intended for healthcare professionals only. For all safety-related information and detailed instructions on the use of this product, please consult the electronic Instructions for Use (eIFU) available at [labeling.merete.de](https://labeling.merete.de). A printed version of the eIFU can be provided free of charge upon request. Ensure that you have read and understood the most current product information before using this product. If you have further questions or in case, you need further information please contact your local Meretes® sales representative or visit our website [merete.de](https://merete.de). Each healthcare professional must use their independent clinical judgment in the diagnosis and treatment of individual patients. Not all patients are suitable for this product or procedure. Check for local product clearances and follow the product-specific and safety-related information of this product.

©2025 Merete GmbH

Follow us on



LinkedIn



Instagram